

# CS STANDARD PRODUCTS


## TECHNICAL CORKS



### CS CUVÉE 0+2


CS Cuvée 0+2 offers superior extraction and compression performance to guarantee the effervescence of the finest Champagne and sparkling wines of the world.

CS Cuvée 0+2 has an agglomerated molded body topped off with two natural cork discs at the end that comes into contact with the wine. Each natural cork disc undergoes the Innocork® process, our patented TCA extraction system, for consistent sensory performance.

	TEST	METHOD	SPECIFICATION	AQL
<b>SENSORY PROPERTIES</b>	<b>RELEASABLE TCA<sup>(1)</sup></b>	ISO 20752	≤ 1.2 ng/L	N.A.
<b>DIMENSIONAL PROPERTIES</b>	<b>LENGTH (STANDARD)</b>	ISO 9727-1; ISO 2859; ISO 4710	<b>48</b> +/- 0.5 mm	<b>4.0</b>
	<b>DIAMETER (STANDARD)</b>	ISO 9727-1; ISO 2859; ISO 4710	<b>29.5 / 30.5 / 31</b> +/- 0.3 mm	<b>1.5</b>
<b>PHYSICAL PROPERTIES</b>	<b>GRANULE SIZE</b>	Internal	<b>3.0 - 7.0 mm</b>	<b>N.A.</b>
	<b>DENSITY</b>	ISO 9727-2; ISO 2859	<b>240 - 300 kg/m<sup>3</sup></b>	<b>4.0</b>
	<b>MOISTURE CONTENT</b>	ISO 9727-3; ISO 2859	<b>&lt; 7.0 %</b>	<b>2.5</b>
	<b>TORSION</b>	NP 2803-6; ISO 2859	<b>≥ 40 daN.cm</b>	<b>N.A.</b>
<b>AVAILABLE WASHES</b>				L04

100% of our natural cork discs undergo the patented INNOCORK® process.

100% of our granules undergo the Vapex® process.

 DS100® service is available.

<sup>(1)</sup> Average value - ISO 20752 | <sup>(2)</sup> Lots > 10,000 corks

N.D. - Not detectable | N.A. - Not applicable | AQL - Acceptable Quality Level

Edition: March 2019